

Snap-On Pet Bandanas



More pictures at www.LoveMendingSouls.org

<u>Snap-on bandana folds over collar, with snaps at each of 2 collar corners. Meant for calming</u> <u>sprays to reduce anxiety and stress behaviors</u>

 Cotton or Poly (see sizing chart). Remember, less material = less heat retention. Measure collar, choose "Across collar" to be half or more collar length.

Approx. collar length: Chihuahua 8"; Whippets 11"; and Golden Retriever 21"

- ♥ 2 sets of plastic snaps (=4 thumb tacks, plus 2 outers and 2 inners)
- ♥ MS flat label and optional fun appliques

Instructions: Use 3/8" seam allowance, back stitch all start and stops.

1. **Collar Wrap:** 6x7.5" Small, 6x10.5" Med and 6x14.5" Large. *Serge all edges, OR attach double bias tape around, OR Hem by folding all edges 1/4" down, press, then fold 1/4" down again, press.

2. Fold lengthways. Attach snaps at four corners, above the hemmed edge. (Thumb tacks on outside, inside fabric has an outer snap to attach to the inner snap.) Test for functionality.

3. **Triangle**: Cut (5x5" Small, 7x7" Med, 10x10" Large) squares diagonally to make two triangles (set aside extra triangle for another bandana). Finish all raw edges as above.* Attach label and appliques on outside by zigzag stitching all around edges to prevent fraying.

4. Find centerpoint of top Triangle edge (can fold in half to find), align to center of one of Collar Wrap's long edge. Sew across, wrong sides together.

Collar Wrap + Triangle = Bandana

Collar Wrap over collar, single layer Triangle attached at center. Showing left set of snaps not closed around collar.

